

C-X-C motif chemokine 13

Catalog No: #AP82206



Package Size: #AP82206-1 50ug #AP82206-2 100ug #AP82206-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-X-C motif chemokine 13
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Angie, B cell-attracting chemokine 1, BCA-1, B lymphocyte chemoattractant, CXC chemokine BLC, Small-inducible cytokine B13, BCA1, BLC, SCYB13
Accession No.	Uniprot ID: O43927
Uniprot	O43927
GeneID	10563;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CXCL13

Note: This product is for in vitro research use only