

# Versican core protein

Catalog No: #AP82207

Package Size: #AP82207-1 50ug #AP82207-2 100ug #AP82207-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Versican core protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Chondroitin sulfate proteoglycan core protein 2,Chondroitin sulfate proteoglycan 2,Large fibroblast proteoglycan,PG-M,CSPG2
Accession No.	Uniprot ID:Q90953
Uniprot	Q90953
GeneID	395565;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:VCAN

Note: This product is for in vitro research use only