

Transcriptional repressor CTCFL

Catalog No: #AP82210



Package Size: #AP82210-1 50ug #AP82210-2 100ug #AP82210-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transcriptional repressor CTCFL
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Brother of the regulator of imprinted sites,CCCTC-binding factor,CTCF paralog,CTCF-like protein,Cancer/testis antigen 27,CT27,Zinc finger protein CTCF-T,BORIS
Accession No.	Uniprot ID:Q8NI51
Uniprot	Q8NI51
GeneID	140690;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CTCFL

Note: This product is for in vitro research use only