

Cytotoxic T-lymphocyte protein 4

Catalog No: #AP82211



Package Size: #AP82211-1 50ug #AP82211-2 100ug #AP82211-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cytotoxic T-lymphocyte protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Cytotoxic T-lymphocyte-associated antigen 4,CTLA-4,CD152
Accession No.	Uniprot ID:Q9XSI1
Uniprot	Q9XSI1
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CTLA4

Note: This product is for in vitro research use only