

H(+)/Cl(-) exchange transporter 5

Catalog No: #AP82214

Package Size: #AP82214-1 50ug #AP82214-2 100ug #AP82214-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	H(+)/Cl(-) exchange transporter 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Chloride channel protein 5, CIC-5, Chloride transporter CIC-5, CLCK2
Accession No.	Uniprot ID:P51795
Uniprot	P51795
GeneID	1184;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CLCN5

Note: This product is for in vitro research use only