

Collectin-11

Catalog No: #AP82215



Package Size: #AP82215-1 50ug #AP82215-2 100ug #AP82215-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Collectin-11
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Collectin kidney protein 1,CL-K1
Accession No.	Uniprot ID:Q17QH6
Uniprot	Q17QH6
GeneID	613991;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:COLEC11

Note: This product is for in vitro research use only