

Clusterin

Catalog No: #AP82216



Package Size: #AP82216-1 50ug #AP82216-2 100ug #AP82216-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Clusterin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	CP40, Complement cytolysis inhibitor, CLI
Accession No.	Uniprot ID:Q29549
Uniprot	Q29549
GeneID	397025;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CLU

Note: This product is for in vitro research use only