

Complement C2

Catalog No: #AP82219



Package Size: #AP82219-1 50ug #AP82219-2 100ug #AP82219-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Complement C2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	C3/C5 convertase
Accession No.	Uniprot ID:Q3SYW2
Uniprot	Q3SYW2
GeneID	515440;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: C2

Note: This product is for in vitro research use only