

Cofilin-1

Catalog No: #AP82221



Package Size: #AP82221-1 50ug #AP82221-2 100ug #AP82221-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cofilin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	18 kDa phosphoprotein,p18,Cofilin, non-muscle isoform,CFL
Accession No.	Uniprot ID:P23528
Uniprot	P23528
GeneID	1072;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CFL1

Note: This product is for in vitro research use only