

Carboxypeptidase N subunit 2

Catalog No: #AP82222

Package Size: #AP82222-1 50ug #AP82222-2 100ug #AP82222-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Carboxypeptidase N subunit 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Carboxypeptidase N 83 kDa chain, Carboxypeptidase N large subunit, Carboxypeptidase N polypeptide 2, Carboxypeptidase N regulatory subunit
Accession No.	Uniprot ID:Q9DBB9
Uniprot	Q9DBB9
GeneID	71756;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cpn2

Note: This product is for in vitro research use only