

# CREB3 regulatory factor

Catalog No: #AP82225



Package Size: #AP82225-1 50ug #AP82225-2 100ug #AP82225-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CREB3 regulatory factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Luman recruitment factor,LRF
Accession No.	Uniprot ID:Q8CDG5
Uniprot	Q8CDG5
GeneID	77128;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Crebrf

Note: This product is for in vitro research use only