

Cystatin-S

Catalog No: #AP82227



Package Size: #AP82227-1 50ug #AP82227-2 100ug #AP82227-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cystatin-S
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cystatin-4,Cystatin-SA-III,Salivary acidic protein 1
Accession No.	Uniprot ID:P01036
Uniprot	P01036
GeneID	1472;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CST4

Note: This product is for in vitro research use only