

Macrophage colony-stimulating factor 1 receptor

Catalog No: #AP82229



Package Size: #AP82229-1 50ug #AP82229-2 100ug #AP82229-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Macrophage colony-stimulating factor 1 receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CSF-1 receptor,CSF-1-R,Proto-oncogene c-Fms,CD115,Csfmr,Fms
Accession No.	Uniprot ID:P09581
Uniprot	P09581
GeneID	12978;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Csf1r

Note: This product is for in vitro research use only