

Casein kinase II subunit alpha

Catalog No: #AP82230



Package Size: #AP82230-1 50ug #AP82230-2 100ug #AP82230-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Casein kinase II subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CK II alpha, Ckii
Accession No.	Uniprot ID:Q60737
Uniprot	Q60737
GeneID	12995;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Csnk2a1

Note: This product is for in vitro research use only