

Choline transporter-like protein 4

Catalog No: #AP82233



Package Size: #AP82233-1 50ug #AP82233-2 100ug #AP82233-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Choline transporter-like protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Solute carrier family 44 member 4,C6orf29,CTL4,NG22,UNQ441/PRO874
Accession No.	Uniprot ID:Q53GD3
Uniprot	Q53GD3
GeneID	80736;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SLC44A4

Note: This product is for in vitro research use only