

Cubilin

Catalog No: #AP82234



Package Size: #AP82234-1 50ug #AP82234-2 100ug #AP82234-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cubilin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	460 kDa receptor,Intestinal intrinsic factor receptor,Intrinsic factor-cobalamin receptor,Intrinsic factor-vitamin B12 receptor,IFCR
Accession No.	Uniprot ID:O60494
Uniprot	O60494
GeneID	8029;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CUBN

Note: This product is for in vitro research use only