

Sorbitol dehydrogenase

Catalog No: #AP82235



Package Size: #AP82235-1 50ug #AP82235-2 100ug #AP82235-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sorbitol dehydrogenase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	L-iditol 2-dehydrogenase
Accession No.	Uniprot ID:Q29318
Uniprot	Q29318
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SORD

Note: This product is for in vitro research use only