

Endoribonuclease Dicer

Catalog No: #AP82236



Package Size: #AP82236-1 50ug #AP82236-2 100ug #AP82236-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Endoribonuclease Dicer
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Helicase with RNase motif, Helicase MOI, DICER, HERNA, KIAA0928
Accession No.	Uniprot ID:Q9UPY3
Uniprot	Q9UPY3
GeneID	23405;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: DICER1

Note: This product is for in vitro research use only