

## Cytochrome P450 1A2

Catalog No: #AP82238



Package Size: #AP82238-1 50ug #AP82238-2 100ug #AP82238-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Cytochrome P450 1A2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CYP1A2,Cholesterol 25-hydroxylase,Cytochrome P450-P2,Cytochrome P450-P3,Cyp1a-2
Accession No.	Uniprot ID:P00186
Uniprot	P00186
GeneID	13077;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Cyp1a2

Note: This product is for in vitro research use only