

DENN domain-containing protein 2D

Catalog No: #AP82240



Package Size: #AP82240-1 50ug #AP82240-2 100ug #AP82240-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	DENN domain-containing protein 2D
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Accession No.	Uniprot ID:Q9H6A0
Uniprot	Q9H6A0
GeneID	79961;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DENND2D

Note: This product is for in vitro research use only