

Carnosine synthase 1

Catalog No: #AP82245



Package Size: #AP82245-1 50ug #AP82245-2 100ug #AP82245-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Carnosine synthase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	ATP-grasp domain-containing protein 1,ATPGD1
Accession No.	Uniprot ID:D3KCC4
Uniprot	D3KCC4
GeneID	100359387;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CARNS1

Note: This product is for in vitro research use only