

## Platelet basic protein

Catalog No: #AP82246



Package Size: #AP82246-1 50ug #AP82246-2 100ug #AP82246-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Platelet basic protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PBP,C-X-C motif chemokine 7,Leukocyte-derived growth factor,LDGF,Macrophage-derived growth factor,MDGF,Small-inducible cytokine B7,CTAP3,CXCL7,SCYB7,TGB1,THBGB1
Accession No.	Uniprot ID:P02775
Uniprot	P02775
GeneID	5473;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PPBP

Note: This product is for in vitro research use only