

C-X-C chemokine receptor type 7

Catalog No: #AP82247



Package Size: #AP82247-1 50ug #AP82247-2 100ug #AP82247-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-X-C chemokine receptor type 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CXC-R7,CXCR-7,Chemokine orphan receptor 1,G-protein coupled receptor 159,G-protein coupled receptor RDC1 homolog,RDC-1,CMKOR1,GPR159,RDC1
Accession No.	Uniprot ID:P25106
Uniprot	P25106
GeneID	57007;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CXCR7

Note: This product is for in vitro research use only