

Cannabinoid receptor 2

Catalog No: #AP82248



Package Size: #AP82248-1 50ug #AP82248-2 100ug #AP82248-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cannabinoid receptor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CB-2,CX5,CB2A,CB2B
Accession No.	Uniprot ID:P34972
Uniprot	P34972
GeneID	1269;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CNR2

Note: This product is for in vitro research use only