

## Cardiotrophin-like cytokine factor 1

Catalog No: #AP82250



Package Size: #AP82250-1 50ug #AP82250-2 100ug #AP82250-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Cardiotrophin-like cytokine factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	B-cell-stimulating factor 3,BSF-3,Novel neurotrophin-1,NNT-1,Bsf3,Nnt1
Accession No.	Uniprot ID:Q9QZM3
Uniprot	Q9QZM3
GeneID	56708;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Clcf1

Note: This product is for in vitro research use only