

## Calcium-binding mitochondrial carrier protein Aralar1

Catalog No: #AP82251



Package Size: #AP82251-1 50ug #AP82251-2 100ug #AP82251-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Calcium-binding mitochondrial carrier protein Aralar1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Mitochondrial aspartate glutamate carrier 1, Solute carrier family 25 member 12, Aralar1
Accession No.	Uniprot ID:Q8BH59
Uniprot	Q8BH59
GeneID	78830;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Slc25a12

Note: This product is for in vitro research use only