

## C-X-C motif chemokine 2

Catalog No: #AP82252



Package Size: #AP82252-1 50ug #AP82252-2 100ug #AP82252-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | C-X-C motif chemokine 2  |
| Brief Description  | Recombinant Protein  |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE   |
| Species Reactivity | Hu   |
| Other Names        | Growth-regulated protein beta,Gro-beta,Macrophage inflammatory protein<br>2-alpha,MIP2-alpha,GRO2,GROB,MIP2A,SCYB2 |
| Accession No.      | Uniprot ID:P19875  |
| Uniprot            | P19875   |
| GeneID             | 2920;  |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)   |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CXCL2

Note: This product is for in vitro research use only