

Cytochrome c oxidase subunit 5A, mitochondrial

Catalog No: #AP82254



Package Size: #AP82254-1 50ug #AP82254-2 100ug #AP82254-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytochrome c oxidase subunit 5A, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Cytochrome c oxidase polypeptide Va
Accession No.	Uniprot ID:P99501
Uniprot	P99501
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:COX5A

Note: This product is for in vitro research use only