

## Dihydropyrimidinase-related protein 2

Catalog No: #AP82258



Package Size: #AP82258-1 50ug #AP82258-2 100ug #AP82258-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Dihydropyrimidinase-related protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	DRP-2,Collapsin response mediator protein 2,CRMP-2,Turned on after division 64 kDa protein,TOAD-64
Accession No.	Uniprot ID:P47942
Uniprot	P47942
GeneID	25416;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Dpysl2

Note: This product is for in vitro research use only