

Protein ripply3

Catalog No: #AP82259



Package Size: #AP82259-1 50ug #AP82259-2 100ug #AP82259-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein ripply3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Down syndrome critical region protein 6,DSCR6
Accession No.	Uniprot ID:P57055
Uniprot	P57055
GeneID	53820;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RIPPLY3

Note: This product is for in vitro research use only