

Desmoglein-3

Catalog No: #AP82260

Package Size: #AP82260-1 50ug #AP82260-2 100ug #AP82260-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Desmoglein-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	130 kDa pemphigus vulgaris antigen homolog
Accession No.	Uniprot ID:O35902
Uniprot	O35902
GeneID	13512;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Dsg3

Note: This product is for in vitro research use only