

## C-X-C motif chemokine 6

Catalog No: #AP82264



Package Size: #AP82264-1 50ug #AP82264-2 100ug #AP82264-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	C-X-C motif chemokine 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Chemokine alpha 3,CKA-3,Granulocyte chemotactic protein 2,GCP-2,Small-inducible cytokine B6,GCP2,SCYB6
Accession No.	Uniprot ID:P80162
Uniprot	P80162
GeneID	6372;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CXCL6

Note: This product is for in vitro research use only