

C-X-C motif chemokine 5

Catalog No: #AP82266

Package Size: #AP82266-1 50ug #AP82266-2 100ug #AP82266-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-X-C motif chemokine 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cytokine LIX, Small-inducible cytokine B5, Scyb5
Accession No.	Uniprot ID:P97885
Uniprot	P97885
GeneID	60665;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cxcl5

Note: This product is for in vitro research use only