

C-X-C chemokine receptor type 2

Catalog No: #AP82267



Package Size: #AP82267-1 50ug #AP82267-2 100ug #AP82267-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-X-C chemokine receptor type 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CXC-R2,GRO/MGSA receptor,High affinity interleukin-8 receptor B,IL-8R B,CD182,Cmkar2,Gpcr16,Il8rb
Accession No.	Uniprot ID:P35343
Uniprot	P35343
GeneID	12765;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cxcr2

Note: This product is for in vitro research use only