

Chorionic somatomammotropin hormone 1

Catalog No: #AP82268



Package Size: #AP82268-1 50ug #AP82268-2 100ug #AP82268-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Chorionic somatomammotropin hormone 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	BPLP-I,Placental lactogen I,PL1
Accession No.	Uniprot ID:P09611
Uniprot	P09611
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CSH1

Note: This product is for in vitro research use only