

Cubilin

Catalog No: #AP82270



Package Size: #AP82270-1 50ug #AP82270-2 100ug #AP82270-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cubilin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Intrinsic factor-cobalamin receptor,Ifcr
Accession No.	Uniprot ID:Q9JLB4
Uniprot	Q9JLB4
GeneID	65969;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cubn

Note: This product is for in vitro research use only