

Dynamin-1-like protein

Catalog No: #AP82271



Package Size: #AP82271-1 50ug #AP82271-2 100ug #AP82271-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Dynamin-1-like protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Dynamin family member proline-rich carboxyl-terminal domain less,Dymple,Dynamin-related protein 1,Drp1
Accession No.	Uniprot ID:Q8K1M6
Uniprot	Q8K1M6
GeneID	74006;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Dnm1l

Note: This product is for in vitro research use only