

Dynamin-1-like protein

Catalog No: #AP82272



Package Size: #AP82272-1 50ug #AP82272-2 100ug #AP82272-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dynamin-1-like protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Dynamin-like protein,Dlp1
Accession No.	Uniprot ID:O35303
Uniprot	O35303
GeneID	114114;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Dnm1l

Note: This product is for in vitro research use only