

Dihydropyrimidinase-related protein 3

Catalog No: #AP82273



Package Size: #AP82273-1 50ug #AP82273-2 100ug #AP82273-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dihydropyrimidinase-related protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	DRP-3,Collapsin response mediator protein 4,CRMP-4,Unc-33-like phosphoprotein 1,ULIP-1,CRMP4,DRP3,ULIP,ULIP1
Accession No.	Uniprot ID:Q14195
Uniprot	Q14195
GeneID	1809;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DPYSL3

Note: This product is for in vitro research use only