

C-X-C motif chemokine 9

Catalog No: #AP82274



Package Size: #AP82274-1 50ug #AP82274-2 100ug #AP82274-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-X-C motif chemokine 9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Gamma-interferon-induced monokine, Monokine induced by interferon-gamma, MIG, Protein m119, Small-inducible cytokine B9, Mig, Scyb9
Accession No.	Uniprot ID:P18340
Uniprot	P18340
GeneID	17329;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cxcl9

Note: This product is for in vitro research use only