

C-X-C chemokine receptor type 2

Catalog No: #AP82275



Package Size: #AP82275-1 50ug #AP82275-2 100ug #AP82275-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-X-C chemokine receptor type 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	CXC-R2,CXCR-2,GRO/MGSA receptor,High affinity interleukin-8 receptor B,IL-8R B,IL8RB,CD182
Accession No.	Uniprot ID:O97571
Uniprot	O97571
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CXCR2

Note: This product is for in vitro research use only