

DnaJ homolog subfamily B member 1

Catalog No: #AP82282



Package Size: #AP82282-1 50ug #AP82282-2 100ug #AP82282-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | DnaJ homolog subfamily B member 1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | DnaJ protein homolog 1,Heat shock 40 kDa protein 1,HSP40,Human DnaJ protein 1,hDj-1,DNAJ1,HDJ1,HSPF1 |
| Accession No. | Uniprot ID:P25685 |
| Uniprot | P25685 |
| GeneID | 3337; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DNAJB1

Note: This product is for in vitro research use only