

# Carnosine synthase 1

Catalog No: #AP82285



Package Size: #AP82285-1 50ug #AP82285-2 100ug #AP82285-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Carnosine synthase 1                                   |
| Brief Description  | Recombinant Protein                                    |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE                           |
| Species Reactivity | Hu   |
| Other Names        | ATP-grasp domain-containing protein 1,ATPGD1,KIAA1394  |
| Accession No.      | Uniprot ID:A5YM72                                      |
| Uniprot            | A5YM72   |
| GeneID             | 57571;   |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.) |

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: CARNS1

Note: This product is for in vitro research use only