

Alpha-crystallin B chain

Catalog No: #AP82288



Package Size: #AP82288-1 50ug #AP82288-2 100ug #AP82288-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Alpha-crystallin B chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Alpha(B)-crystallin
Accession No.	Uniprot ID:P23928
Uniprot	P23928
GeneID	25420;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cryab

Note: This product is for in vitro research use only