

Cytochrome b5

Catalog No: #AP82295



Package Size: #AP82295-1 50ug #AP82295-2 100ug #AP82295-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytochrome b5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Microsomal cytochrome b5 type A, MCB5, CYB5
Accession No.	Uniprot ID:P00167
Uniprot	P00167
GeneID	1528;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CYB5A

Note: This product is for in vitro research use only