

Cartilage acidic protein 1

Catalog No: #AP82297



Package Size: #AP82297-1 50ug #AP82297-2 100ug #AP82297-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cartilage acidic protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	68 kDa chondrocyte-expressed protein, CEP-68, ASPIC, ASPIC1, CEP68
Accession No.	Uniprot ID:Q9NQ79
Uniprot	Q9NQ79
GeneID	55118;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CRTAC1

Note: This product is for in vitro research use only