

Dihydropteridine reductase

Catalog No: #AP82300



Package Size: #AP82300-1 50ug #AP82300-2 100ug #AP82300-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dihydropteridine reductase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	HDHPR, Quinoid dihydropteridine reductase, Short chain dehydrogenase/reductase family 33C member 1, DHPR, SDR33C1
Accession No.	Uniprot ID:P09417
Uniprot	P09417
GeneID	5860;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0) gene_full_name: QDPR

Note: This product is for in vitro research use only