

D(2) dopamine receptor

Catalog No: #AP82304

Package Size: #AP82304-1 50ug #AP82304-2 100ug #AP82304-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	D(2) dopamine receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Dopamine D2 receptor
Accession No.	Uniprot ID:Q9GJU1
Uniprot	Q9GJU1
GeneID	403701;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: DRD2

Note: This product is for in vitro research use only