

Dynamin-1

Catalog No: #AP82305



Package Size: #AP82305-1 50ug #AP82305-2 100ug #AP82305-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dynamin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Dnm,Kiaa4093
Accession No.	Uniprot ID:P39053
Uniprot	P39053
GeneID	13429;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Dnm1

Note: This product is for in vitro research use only