

C-X-C motif chemokine 10

Catalog No: #AP82306



Package Size: #AP82306-1 50ug #AP82306-2 100ug #AP82306-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description

|                    |                                                                                                                                                      |
|--------------------|------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name       | C-X-C motif chemokine 10                                                                                                                             |
| Brief Description  | Recombinant Protein                                                                                                                                  |
| Host Species       | E.coli                                                                                                                                               |
| Purification       | Greater than 90% by SDS-PAGE                                                                                                                         |
| Species Reactivity | Ms                                                                                                                                                   |
| Other Names        | 10 kDa interferon gamma-induced protein,Gamma-IP10,C7,Interferon-gamma induced protein<br>CRG-2,Small-inducible cytokine B10,Crg2,Ifi10,Inp10,Scyb10 |
| Accession No.      | Uniprot ID:P17515                                                                                                                                    |
| Uniprot            | P17515                                                                                                                                               |
| GeneID             | 15945;                                                                                                                                               |
| Tag Info           | His                                                                                                                                                  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)                                                                                               |

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Cxcl10

Note: This product is for in vitro research use only